

Spinal cord stimulation (HF-SCS) at 10 kHz for the treatment of chronic focal neuropathic post-surgical pain

Mayank Gupta MD¹; James Scowcroft MD²; Daniel Kloster MD²; Maged Guirguis MD³; Jonathan Carlson MD⁴; Tory McJunkin MD⁴; Gassan Chaiban MD³; Atef Israel MD¹; Jeyakumar Subbaroyan PhD⁵
¹Menorah Medical Center, Overland Park, MO; ²KC Pain Centers, Lee's Summit, MO; ³Ochsner Health System, New Orleans, LA; ⁴AZ Pain Specialists, Scottsdale, AZ; ⁵Nevro Corp., Redwood City, CA

Introduction

Chronic post-surgical pain (CPSP) is one of the largest sub-groups of focal mononeuropathies. Incidence of chronic pain from iatrogenic peripheral nerve injuries resulting from surgeries like mastectomy, thoracotomy and amputation vary from 5-85%¹. The focal nature of the pain may warrant targeting a neural structure like dorsal root ganglion (DRG) to elicit pain relief. However, early pre-clinical evidence from high frequency SCS (HF-SCS) at 10 kHz demonstrated an inhibitory effect on projection interneurons in the superficial layers of dorsal horn. Thus, we hypothesized that HF-SCS at 10 kHz may provide effective pain relief in focal CPSP conditions.

Methods

- Major inclusion criteria
 - Subjects with chronic, intractable pain of ≥ 5 cm (on a 0-10 cm visual analog scale [VAS]) of the trunk, upper or the lower limb from CPSP
 - A score of ≥ 4 in the Douleur Neuropathique 4 (DN4) questionnaire
- Major exclusion criteria: Significant spinal stenosis, epidural scarring or symptoms of myelopathy
- Subjects enrolled following Institutional Review Board approval
- Each subject implanted with two epidural leads spanning appropriate vertebral bodies as determined by the location of pain
- Senza system (Nevro Corp., Redwood City, CA) implanted in subjects with successful trial stimulation ($\geq 40\%$ pain relief)
- Safety and effectiveness endpoints captured up to 12 months post-implant
- Interim three month results presented (mean \pm standard error) in the permanent implant population

Figure 1 (Left). Lead location for a subject with neuropathic groin pain post-hernia repair. The subject's baseline pain was 8.2 cm on VAS. The lead placement was anatomically guided and there was no intra-operative paresthesia mapping involved.






Figure 2. Representative pain maps in a subject with pain post-hernia repair (middle) and pain post-knee surgery (right).


Results: Trial Success

- Trialed: 28
- Trial success: 23/28 (82%)
- Location of pain
 - Upper extremity: 3
 - Trunk: 8
 - Lower extremity: 17

Results: Pain Scores




Results: Pain Disability Index (PDI)



- Disability assessed by PDI (Range: 0-70)
- Lower index score = Less disability
- Minimal clinically important difference (MCID): 8.5-9.5²
- At 3 months: 20.7 point reduction in PDI or over 2X MCID

Results: Pain Interference



- Measured by McGill Pain Questionnaire (SF-MPQ-2)
- Significant reduction in all dimensions of pain including affective descriptors

References

- Macrae WA (2008). Chronic post-surgical pain. Br J Anaesth. 101(1): 77-86.
- Soer et al. (2012) Spine.